


Option 1: MedImage Server on each PC (AtomJump's proxy server used)

1. Phone snaps photo, with name
2. AtomJump cloud receives uploaded photo
3. User's PC downloads photo
4. User uploads named photo onto application



Pros:

- No software to be written.
- Simple setup - only a standard Windows installer need be installed on each user's PC.

Cons:

- Installation on each PC must likely be agreed, and potentially, completed by IT
- Only 5GB of data transfer is free, after which time there is an ongoing data transfer cost.
- Double step phone-to-PC, PC-to-application user process