


### Option 3: MedImage Proxy Server on application server.

1. Phone snaps photo, with name
2. Application server, with MedImage proxy software server receives uploaded photo
3. 'Glue' code transfers named photo directly into application



#### Pros:

- Direct phone-to-system upload. No additional upload step is required on the user's desktop.
- IT does not need to be involved (only application developers)

#### Cons:

- Depending on the image upload file-system on the application, some 'glue' code may need to be written by the application developers to get the images into the correct format/directories. (i.e. some initial setup cost)